W10

Instrucciones / Instructions / Instruções / Instructions /
Istruzioni / Anweisungen / Instructies / Návod / Pokyny

ES

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Keep your purchase receipt in order to exercise your consumer rights. The company reserves the right to make modifications to the product with respect to what is described in this instruction manual. Malfunctions caused by natural disasters, maintenance carried out by a non-original factory, or unauthorized modifications are not covered by the warranty.

Designed to assist people with functional mobility limitations who require safe support while walking.

For people with visual impairments, dyslexia, or cognitive disabilities, please share the instructions with family members.

Before use, make sure all functions operate properly. Stop using the walker if it is damaged or malfunctioning and contact the authorized distributor. This walker is intended for indoor use and for users with a maximum weight of 100 kg.

The maximum load weight in the basket is 5 kg.

Be careful and use the braking system when walking on ramps.

Do not store in direct sunlight or in places with heat/humidity to avoid malfunction or deformation.

Do not use on ramps or uneven surfaces on rainy days or slippery floors. Improper use may cause dangerous situations: do not use the product to transport a person.

Please refer to the Warranty Terms and Conditions at www.playcare.es

Washing instructions:

Read the washing label carefully.

Machine wash on a warm gentle cycle. Line dry.

Do not use solvents or bleach that may damage the materials.

Metal and plastic parts may be hand washed with warm water and mild soap. Do not use aggressive substances.

Dry the walker frame carefully after contact with water (or rain) to prevent moisture.

Do not modify or disassemble the walker frame to avoid damage or malfunction.

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EN

Once the walker has reached the end of its useful life, it must be disposed of in accordance with the environmental regulations applicable in your municipality or autonomous community. To facilitate recycling, it is recommended to disassemble the walker by separating its components (metal frame, wheels, handles, etc.). You can contact the nearest recycling point or urban waste facility for information on how to properly manage these materials. Packaging materials (cardboard, plastics, protections) should also be placed in the appropriate recycling containers or returned to the distributor for proper handling.

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(01) 04713361521683

Class I

CHARMING WORLD INC.
3F-1, No.75, Sec. 3, Blade Road, Songshan
District, Taipei City, 10554, Taiwan

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KOMAS Medical Technology GmbH,
PLAY S.A.